


# GELLI-ROLL®

*Reusable Warming and Cooling with Comfort*



**THE GELLI-ROLL COMBINES THE ABILITY TO WARM AND COOL YOUR PATIENT WITH COMFORT.**

- Allows complete access to the patient
- Easy to clean with disinfectants
- Reduces set-up and breakdown time
- Supports most patients, infant to bariatric

The revolutionary and simple-to-use Gelli-Roll® combines patient temperature management and comfort. The Gelli-Roll® is a reusable water blanket encapsulated in Akton Polymer. It may be used before, during, or after surgery to effectively keep your patient normothermic.

Many medical professionals are looking for ways to cut costs and eliminate environmental impact. In most cases, it pays for itself in a matter of months.


Delivering Tomorrow's Temperature  
Management Solutions Today

**CSZ**  
Cincinnati Sub-Zero

# GELLI-ROLL®

*Reusable Warming and Cooling with Comfort*

Ordering Information		
Cat. No.	Description	Quantity
195P	Gelli-Roll® - Adult 73.5" x 21.5" x 0.6" (186.7cm x 54.6cm x 1.6cm)	ea.
194P	Gelli-Roll® - Pediatric 31.3" x 23.7" x 0.6" (79.4cm x 60.3cm x 1.6cm)	ea.
193P	Gelli-Roll® - Infant 20" x 13" x 0.6" (50.8cm x 33cm x 1.6cm)	ea.
Hyper-Hypothermia units		
233	Blanketrol® III	115v, 230v
222R	Blanketrol® II (not available in Europe)	100v, 115v, 230v, 240v
111W	Norm-O-Temp®	100v, 115v, 230v, 240v
Reusable Connecting Hose (must purchase 286 hose to use Gelli-Roll with hyper-hypothermia units)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.



## Product Specifications:

### Dimensions:

#### Adult:

73.5in x 21.5in x 0.6in  
(186.7cm x 54.61cm x 1.6cm)

**Weight:** 35 lb. (15.9 kg)

#### Pediatric:

31.3in x 23.7in x 0.6in  
(79.38cm x 60.3cm x 1.6cm)

**Weight:** 25 lb (11.34 kg)

#### Infant:

20in x 13in x 0.6in  
(50.8cm x 33.02cm x 1.6cm)

**Weight:** 15 lb (6.8 kg)

**CSZ**  
Cincinnati Sub-Zero

Phone: 513-772-8810  
Toll Free: 800-989-7373  
Fax: 513-772-9119

[www.cszmedical.com](http://www.cszmedical.com)

© 2005 by Cincinnati Sub-Zero Products, Inc.  
99202 6/09