

KOOL-KIT[®]

Therapeutic Temperature Management System



WHEN IT'S TIME TO COOL, CSZ OFFERS AN EFFECTIVE SOLUTION WITH IT'S KOOL-KIT[®].

- About 7% of blood flows through the area treated by the head wrap every minute.¹
- About 33% of blood flows through the vasculature exclusively treated by the Kool-Kit[®] vest.¹
- About 26% of blood flows through the groin area treated by the lower body blanket.¹
- Can cover approximately 66% of your patient, allowing for excellent body surface coverage.¹

With the Kool-Kit[®] caregivers are able to have easy access to their patients, maximum surface-area coverage, and finely administer therapeutic cooling and re-warming for each patient. The CSZ Kool-Kit[®] combines a head wrap, patient vest, and lower body blanket with the advanced Blanketrol[®] III hyper-hypothermia system for an integrated, modular approach to therapeutic temperature management.

1. Data on file at Cincinnati Sub-Zero

Cincinnati Sub-Zero has been managing your patient temperature needs for over 45 years.

CSZ
Cincinnati Sub-Zero

KOOL-KIT®

Therapeutic Temperature Management System

Ordering Information		
Cat. No.	Description	Quantity
900	Kool-Kit® (1 head wrap, 1 patient vest & 1 lower body blanket)	5/box
910	Kool-Kit® Head Wrap (1 head wrap & 1 lower body blanket)	5/box
920	Kool-Kit® Patient Vest (1 patient vest & 1 lower body blanket)	5/box
950	Kool-Kit® Neonate (1 maxi-therm lite blanket, 1 9Fr general purpose probe, 1 set pink mittens, 1 set blue mittens, 1 set white socks) (NOT AVAILABLE INTERNATIONALLY)	5/box
Connecting Hoses (Must purchase one, two or three hoses depending on the kit that is purchased)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.
Hyper-Hypothermia units		
233	Blanketrol® III	115v or 230v

Key Benefits

- Single Patient-Use
- Allows for easy access to patient & monitoring devices
- Cost effective
- Larger body surface area coverage with no adhesive



Product Specifications:

Head Wrap Dimensions:
16.4in x 24.6in (41.7cm x 62.5cm)

Patient Vest Dimensions:
53.2in x 24.9in (135.1cm x 63.2cm)

Lower Body Blanket Dimensions:
33in x 25in (83.8cm x 63.5cm)



CSZ
Cincinnati Sub-Zero

Phone: 513-772-8810
Toll Free: 800-989-7373
Fax: 513-772-9119

www.cszmedical.com

© 2009 by Cincinnati Sub-Zero Products, Inc.
99196 9/10